

# The Hill Cervical Flexion Headpiece

*with Lateral Motion*



To begin cervical flexion, first secure the straps around the occipital region.



Disengage the tilting lock so that the headpiece can flex freely up and down.



Disengage the lateral lock to allow for gentle circumduction and flexion.



Only a 1/2" to 1" of flexion is needed to create a good amount of traction.

For full product details, pricing and a demonstration video visit:  
<http://www.hilllabs.com/chiropractic/Hill-Air-Flex-Table.php>

  
**Hill Laboratories**  
COMPANY

P.O. Box 2028, 3 Bacton Hill Road, Frazer, PA 19355  
1-877-445-5020 • 610-644-2867 • fax 610-647-6297 • [www.hilllabs.com](http://www.hilllabs.com)