

Hill Laboratories Medical Chair

The Versatile and Affordable Chair
for Health Care Professionals



- Shown with Gravel upholstery
and a Beige base.

Versatility

Hill Medical Chair

Working closely with doctors, Hill Laboratories has designed the ideal chair for your practice. From the standard model to one featuring many options, the Hill Medical Chair provides versatility at an affordable price. We will also customize the chair to meet your specific needs.



- Standard Model shown with Indigo upholstery and Beige base.

Adjustable Foot Section Lever



Paper retainer on headpiece.



Paper roll holder on foot cushion.



Standard ★ Features

- 22"-35" high, 24" or 26" wide cushions
- Electric lift-back with radius top
- Flat cushions - for patient comfort, 1" soft foam overlaying 1 1/2" medium density foam
- Full-width foot cushion with manual adjustable tilt
- Paper retainer on headpiece
- Paper roll holder on foot cushion
- Choice of 20 vinyl colors
- Beige or gray powder-coated base with thermo-plastic skirting
- One-year warranty on parts & labor

Custom Options

For ultimate patient comfort and positioning, the Hill Medical Chair can be equipped with an electric foot section and electric Tilt (Trendelenburg) feature. The optional One-Touch Foot Pedal makes it easy to automatically raise or lower the chair with a simple and single “tap”.

Electric Options



75

Electric-Tilt shown with Charcoal upholstery and Gray base.



76A
*Removable
Fold-down Armrests*

Additional



Examining
Shown with full range of motion

Below: Fully reclined with armrests removed



Matching stool (mauve) with or without back



61B
Tapered Back with Adjustable Tilting Headrest



91D
Adjustable Headrest with Face Cradle

For prone positioning

Cushion Options



74

Electric-Foot shown with Teal upholstery and Beige base.



61B

Tapered-Back with headpiece shown with Indigo upholstery and Gray base.



61A

Contour-Seat shown with Suede upholstery and Beige base.



al ★ Options

64

Rotation Armboard

motion from side (left) and end (right)



66

Gooseneck Halogen Lamp

Mounts anywhere on table



**Headrest
Cut-Out**



92H


Small Adjustable Head Support



65 **Receptacle in Base**



Specifications



Powered by an electric actuator 60 Hz, 110 volts

Amperage:

One Motor Table - 2.4, Two Motor Table - 5.1
Three Motor Table - 7.5, Four Motor Table - 10.2

Lifting Capacity:

400 lbs. (600 lbs. capacity upon request)

Height Range:

90ALB - 23" to 36" or 21" to 31"
90D - 22 to 35. With Trendelenburg - 24" to 37" or optional 22" to 32"

Table Base:

Beige or Gray telescoping base.
All metal parts are powder coated.

Shipping Weight:

90ALB - 225 lbs., 90D - 300 lbs.

Warranty: One year

Color ★ Choices

Standard Colors: Cranberry, Oxblood, Plum, Shadow, Mauve, Rose, Crocus, Almond, Suede, Spice, Chestnut, Slate, Hunter, Seafoam, Teal, Blue, Navy, Azure, Gray, Charcoal.

Base Color - Gray: Cranberry, Oxblood, Plum, Shadow, Mauve, Rose, Crocus, Almond, Suede, Spice, Chestnut, Slate, Hunter, Seafoam, Teal, Blue, Navy, Azure, Gray, Charcoal.

Base Color - Black: Cranberry, Oxblood, Plum, Shadow, Mauve, Rose, Crocus, Almond, Suede, Spice, Chestnut, Slate, Hunter, Seafoam, Teal, Blue, Navy, Azure, Gray, Charcoal.

Base Color - Beige: Cranberry, Oxblood, Plum, Shadow, Mauve, Rose, Crocus, Almond, Suede, Spice, Chestnut, Slate, Hunter, Seafoam, Teal, Blue, Navy, Azure, Gray, Charcoal.

Upgraded Finesse™ Upholstery Colors: Bordeaux, Gravel, Iceberg, Hershey, Pine, Colonial Blue, Indigo, Cinder.

Upholster your chair in durable Boltaflex® with Prefix™, scientifically formulated to give superior protection from staining and abrasions (top row). Or, for more of an upscale look and feel, choose puncture-proof Finesse™ material (bottom row). Custom colors are available at additional cost. (Shades may vary from actual material. Colors may change without notice.)

Since 1945 
Hill Laboratories
COMPANY

P.O. Box 2028, 3 N. Bacton Hill Road, Frazer, PA 19355 • Toll-Free 1-877-445-5020
610-644-2867 • fax 610-647-6297 • www.HillLabs.com