

# HillAdjustable HA90MD



*Standard Model  
HA90MD*

*Shown with optional  
contour upholstery,  
adjustable headrest and  
fold-down armrests*

## **HA90MD**

The HA90MD has a host of features designed for the medical profession in addition to 22" to 35" variable height.

Powered by an electric foot pedal, the lift-back section smoothly converts from a full-recline to chair-examination posture. Like the HA90E, the 90MD is equipped with fully adjustable, retractable foot stirrups, 2 1/2" of foam and 24", 26" or optional 30" widths.

Model 90MD has a pull-out, stainless steel tray that slides neatly under the seat when not in use and may be removed for easy cleaning. An optional adjustable-angle seat cushion can be also added providing pelvic tilt. The foot section manually adjusts to any angle for efficient and comfortable positioning. Shorter legs on the base permit easy movement around the table.

Like all Hill Adjustables, mechanical parts are safely enclosed in a grey or beige thermo-plastic base.

# 90MD Options

## Trendelenburg

The 90MD is available with an electric tilt feature for use with Trendelenburg treatments. (pull-out tray not available with this option).



## Exam Stools

Contemporary design with gas-powered height, controlled with a lever, from 19 1/2" - 24 1/2". Available with or without back and upholstered to match your Hill exam table or chair.



## Electric Foot

Stops at any angle and controlled by a foot pedal (pull-out tray not available with this option).



## Armboard

The examination armboard has receiving brackets on both sides of the table.



## Angled Seat

Position your patient with the tilting seat that adjusts manually to two angles.



## Halogen Lamp

With a 34" arm and removable stem which attaches to the table top (includes single electric receptacle).

