

The Counterstrain Table

by Hill Laboratories





from \$4190

Featuring-


- Tilting Headpiece - 90° negative and 30° positive
- Electric Tilt-Back
- Slide-out pelvic section
- Lift leg section - gas spring assisted
- Drop-foot section, 8", Mechlok® operated, 90° negative and 30° positive.
(If omitted, deduct \$125.00)
- Double locks (Mechlok) on each side of the table to raise and lower head and leg section
- Two foot pedals mounted near the front of the table to control head/thoracic area
- Two foot pedal mounted near the base of the table to control height
- Headpiece open at forehead area
- Armrests mounted on head section
- Electric hand-operated switch to elevate table at head end
- 78" long, 24" wide
- 20 designer colors available for the top and gray, beige and black for the base


Hill Laboratories
COMPANY


Hill Counterstrain Osteopathic Table




Motorized Elevation 22" - 35"




Motorized Thoracic Flexion to 85°




Headpiece Adjusts through 135°



Motorized Elevation 22"-35"



Leg Section Flexes to 75°



Wide Range of Positions