

THE HILL ADJUSTABLE MODEL HA90V

featuring electric tilt, electric adjustable height and manual lift-back



For:

- Endovenous Treatment
- Vascular Diagnostic Testing
- Echocardiogram Testing

The vascular exam table Model HA90V is a versatile table offering a sturdy welded-steel base construction and smooth electric elevation.

Standard Features:

- Electric Adjustable Height Base (range 24" to 34" or 26" to 39") controlled by foot pedal
- Electric Tilting Table Top - (head and foot end tilt 20°) also controlled by foot pedal
- Manual Lift Back with switch-activated gas-spring
- Adjustable Leveling Feet
- Top thickness - 2 1/2" medium density foam . Top size - 27" x 63"
- Lifting Capacity - 500 lbs. • Approximate weight - 225 lbs.
- Choice of 20 upholstery colors with gray or beige base

Optional Features:

- | | |
|-------------------------------------------------------|------------------------------------------------------------------------|
| 94V Electric Lift-Back | 49 Foot Board |
| 98 30" Extra Width | 50 Protractor |
| 88 32" Extra Width | 140 Upgrade Upholstery,
Puncture-Resistant Tricot
Finesse |
| 91 10" Face Cut-out | 93 Caster System for mobility.
Engaged with foot lever |
| 91A Insert for Face Cut-out | 96 Paper Roll Bracket |
| 60 Fold-down Echo Window
(* see note below) | |
| 48 Seat Belt | |

* The echo-window features an easy-access 8"x11" fold-down section designed for left lateral position echo cardiology (only available with 30" or 32" width)

O P T I O N S



48



94V



50




60

90V with
Echo Window



49



93



Hill Laboratories
COMPANY