

Hill A-F-T

Automatic **Flexion Table**

Automatic Flexion priced
WELL BELOW
competitive models!



*Shown with optional
drops, navy upholstery
and beige base.*

Features, Options and Specifications

- Variable speed motorized flexion
- Electric adjustable height with dual rocker foot pedals
- Tilting Headpiece
- Thoracic abdominal release locks at any angle
- Slide-out ankle support with ankle harnesses
- Pelvic lateral motion with gear lock
- Patient "Gripper Bar" - patient can adjust to most comfortable length
- Patient Shut-Off Switch
- Optional Drops - cervical, thoracic and pelvic

Specifications

- Electric requirement 60 Hz, 110 volts, 4.8 amps. (220-50 cycle also available)
- Lifting capacity - 450 lbs.
- Length - 5' 5" (extends to 6' 2"), with cervical drop - 5' 7" to 6' 4"
- Width - 24"
- Range of height - 21.5" - 29" or 22.5" - 30"
- Medium density foam top. Firm and soft density upon request.
- Choice of 20 decorator colors for upholstered top. Beige or gray telescoping base.
- One-year warranty on all parts and labor
- Shipping weight 300 lbs.

– See other side for more details –

Features and Options of the Hill A.F.T.



Automatic Flexion

Variable speed motorized flexion with convenient speed dial. Speeds from 1-28 rpm.



Slide-Out Ankle Support

Extends 9" and controlled by a lock lever switch. Includes ankle harnesses



Adjustable Height

Controlled with electric foot pedal. Choice of height range - 21.5" to 29" or 22.5" to 30"



Patient Gripper Bar

Patient squeezes lever and slides bar to adjust to comfortable length.



Tilting Headpiece

Hand lever releases headpiece to lock at any angle from 30° positive to 30° negative.



Patient Shut-Off

Equipped with a patient shut-off switch immediately cutting all table power when pressed.



Thoracic Release

Thoracic Abdominal Release locks at any angle down to 5" and up (above horizontal) to 3"



Optional Drops

Optional cervical, thoracic and pelvic drops may be added to the A.F.T. at additional cost.

Optional Headpieces

With each headpiece below, the adjustable patient gripper bar is individually mounted under each armrest.



Flexion Headpiece

For manual flexion. Lateral and long-axis movement. Includes raised headpiece.



CBP Headpiece

Raised and tilt headpiece with translation.



Raised Headpiece

Raises 5" and equipped with tilting headpiece



Flexion Control Panel

The Flexion Control Panel is conveniently located at the end of the flexion cushion and allows you to set the flexion speed from 1 to 28 rpm. A selector switch permits the choice of stopping the table in a horizontal position or at any angled position.


Hill Laboratories
COMPANY

3 N. Bacton Hill Road, Frazer, PA 19355 • P.O. Box 2028, Frazer, PA 19355
Toll-Free 1-877-445-5020 • 610-644-2867 • fax 610-647-6297 • www.hilllabs.com