

# Hill A-F-T

## Automatic **Flexion Table**

Automatic Flexion priced  
**WELL BELOW**  
competitive models!



*Shown with optional  
drops, navy upholstery  
and beige base.*

## Features, Options and Specifications

- Variable speed motorized flexion
- Electric adjustable height with dual rocker foot pedals
- Tilting Headpiece
- Thoracic abdominal release locks at any angle
- Slide-out ankle support with ankle harnesses
- Pelvic lateral motion with gear lock
- Patient "Gripper Bar" - patient can adjust to most comfortable length
- Patient Shut-Off Switch
- Optional Drops - cervical, thoracic and pelvic

### Specifications

- Electric requirement 60 Hz, 110 volts, 4.8 amps. (220-50 cycle also available)
- Lifting capacity - 450 lbs.
- Length - 5' 5" (extends to 6' 2"), with cervical drop - 5' 7" to 6' 4"
- Width - 24"
- Range of height - 21.5" - 29" or 22.5" - 30"
- Medium density foam top. Firm and soft density upon request.
- Choice of 20 decorator colors for upholstered top. Beige or gray telescoping base.
- One-year warranty on all parts and labor
- Shipping weight 300 lbs.

– See other side for more details –

# Features and Options of the Hill A.F.T.



## Automatic Flexion

Variable speed motorized flexion with convenient speed dial. Speeds from 1-28 rpm.



## Slide-Out Ankle Support

Extends 9", and controlled by a lock lever switch. Includes ankle harnesses



## Adjustable Height

Controlled with electric foot pedal. Choice of height range - 21.5" to 29" or 22.5" to 30"



## Patient Gripper Bar

Patient squeezes lever and slides bar to adjust to comfortable length.



## Tilting Headpiece

Hand lever releases headpiece to lock at any angle from 30° positive to 30° negative.



## Patient Shut-Off

Equipped with a patient shut-off switch immediately cutting all table power when pressed.



## Thoracic Release

Thoracic Abdominal Release locks at any angle down to 5" and up (above horizontal) to 3"



## Optional Drops

Optional cervical, thoracic and pelvic drops may be added to the A.F.T. at additional cost.

## Optional Headpieces

With each headpiece below, the adjustable patient gripper bar is individually mounted under each armrest.



### Flexion Headpiece

For manual flexion. Lateral and long-axis movement. Includes raised headpiece.



### CBP Headpiece

Raised and tilt headpiece with translation.



### Raised Headpiece

Raises 5" and equipped with tilting headpiece



## Flexion Control Panel

The Flexion Control Panel is conveniently located at the end of the flexion cushion and allows you to set the flexion speed from 1 to 28 rpm. A selector switch permits the choice of stopping the table in a horizontal position or at any angled position.

  
**Hill Laboratories**  
COMPANY

3 N. Bacton Hill Road, Frazer, PA 19355 • P.O. Box 2028, Frazer, PA 19355  
Toll-Free 1-877-445-5020 • 610-644-2867 • fax 610-647-6297 • www.hilllabs.com