

# Hill Laboratories

## HA90D

The Versatile and Affordable Chair  
for Dermatology



- Shown with Gravel upholstery  
and a Beige base.

# For Dermatology

## HA90D

Working closely with Dermatologists, Hill Laboratories has designed the ideal chair for dermatology. From the standard 90D to one featuring many options, the Hill Dermatology Chair provides versatility at an affordable price. We will also customize the 90D to meet your specific needs.



- Standard Model shown with Indigo upholstery and Beige base.

Adjustable Foot Section Lever



Paper retainer on headpiece.



Paper roll holder on foot cushion.



## Standard ★ Features

- 22"-35" high, 24" or 26" wide cushions
- Electric lift-back with radius top
- Flat cushions - for patient comfort, 1" soft foam overlaying 1 1/2" medium density foam
- Full-width foot cushion with manual adjustable tilt
- Paper retainer on headpiece
- Paper roll holder on foot cushion
- Choice of 20 vinyl colors
- Beige or gray powder-coated base with thermo-plastic skirting
- One-year warranty on parts & labor

# HA90D Options

For ultimate patient comfort and positioning, the HA90D can be equipped with an electric foot section and electric Trendelenburg feature. The optional One-Touch Foot Pedal makes it easy to automatically raise or lower the table with a simple “tap”.

## Electric Options



75

*Electric-Tilt shown with Charcoal upholstery and Gray base.*



76A

**Removable  
Fold-down Armrests**

## Additional



**Examining**

Shown with full range of motion

**Below:** Fully reclined with armrests removed



*Matching stool (mauve) with or without back*

61B

**Tapered Back with Adjustable Tilting Headrest**



91D

**Adjustable Headrest with Face Cushion**

For prone positioning

# Cushion Options



74

Electric-Foot shown with Teal upholstery and Beige base.



61B

Tapered-Back with headpiece shown with Indigo upholstery and Gray base.



61A

Contour-Seat shown with Suede upholstery and Beige base.




# ★ Options

64

Rotation Armboard

Option from side (left) and end (right)



66

Gooseneck Halogen Lamp

Mounts anywhere on table



Headrest  
t-Out



92H

Small Adjustable Head Support



65 Receptacle in Base



# For Aesthetics

## HA90ALB

The 90ALB (Aesthetic Lift Back) is an extremely affordable chair for skin care with 2" of soft foam over 1 1/2" of medium-density for maximum patient comfort. Many standard and optional features make the 90ALB ideal for aesthetic work.

- Shown with  
Bordeaux upholstery  
and Beige base.


### Easy up. Easy down.

The 90ALB features electric elevation from 23" to 36" or 21" to 31". The standard model has an assisted, gas spring back for easy positioning, an adjustable headrest, fold-down removable armrests, contour back and durable, upgraded designer upholstery.

The HA90ALB can fully recline and is also available with an optional electric lift back and a variety of head cushion options.



# Specifications



Powered by an electric actuator 60 Hz, 110 volts

**Amperage:**

One Motor Table - 2.4, Two Motor Table - 5.1  
Three Motor Table - 7.5, Four Motor Table - 10.2

**Lifting Capacity:**

400 lbs. (600 lbs. capacity upon request)

**Height Range:**

90ALB - 23" to 36" or 21" to 31"  
90D - 22 to 35. With Trendelenburg - 24" to 37"  
or optional 22" to 32"

**Table Base:**


Beige or Gray telescoping base.  
All metal parts are powder coated.

**Shipping Weight:**

90ALB - 225 lbs., 90D - 300 lbs.

**Warranty:** One year

## Color Choices



Upholster your table in durable Boltaflex® with Prefix™, scientifically formulated to give superior protection from staining and abrasions (top row). Or, for more of an upscale look and feel, choose puncture-proof Finesse™ material (bottom row). Custom colors are available at additional cost. (Shades may vary from actual material. Colors may change without notice.)

Since 1945   
**Hill Laboratories**  
COMPANY

Cover photo taken at  
Family Dermatology,  
Warrington, PA