



Hands-Free **Ultrasound** *with Muscle Stimulation*



*Saves Hours of Time Compared
to Manual Ultrasound!*

Quick Set-Up. *Time-Tested Results.*

The Hands-Free Ultrasound is designed with 3 crystals and administers pulsating ultrasound at 1 megahertz. The crystals are designed to harmonically vibrate the entire soundhead which covers a treatment area of 3 5/8" diameter (65cm²). Studies have shown that ultrasound with low intensity and longer treatment times is as effective as conventional, hand-held methods. The Hands-Free Ultrasound allows for extended treatment times without the worry of operator error and fatigue associated with conventional ultrasound.

Benefits

- Saves money, personnel time and energy
- Stationary Soundhead allows for safe, controlled, repeatable treatments
- Combination of ultrasound and muscle stimulation directly through the soundhead

Features

One of the most revolutionary features of the HF27 is the reusable *Gel Pad*.

The Gel Pad eliminates the need for messy gels and contours to curved areas of the body.

Ultrasound and Muscle Stimulation will both conduct through the Gel Pad.



Shown with
Gel Pad &
Retaining Ring

Other Features include -

- Accurate duplication of treatments and intensities
- High voltage muscle stimulator: Frequency 1-100, Intensity 0-100
- Selector Switch for muscle stimulation, ultrasound, or combination through the soundhead
- Safety Start/Stop Button that will not allow operation of any functions until all switches are reset.
- Patient Cut-Off Switch
- Complete with cart and storage tray, adjustable arm, gel pad, spray bottle, weight bag and 2 adjustable straps, and one set of carbon electrodes

See other side for easy set-ups of common treatment areas

Hands-Free Applications of the HF27

Pictured below are several common treatment setups for the HF27; each takes less than a minute to prepare.



Wrist



Elbow



Shoulder



Upper Back



Neck



Lower Back



Bottom of Foot



Side of Foot



Knee



HF27
Soundhead
Setup

Soundhead
Retainer Ring
Gel Pad
Wet Paper Towel

Distributed by:


Hill Laboratories
COMPANY

3 N. Bacton Hill Road, Frazer, PA 19355

Toll-Free 1-877-445-5020

Ph. 610-644-2867 • Fax: 610-647-6297

www.hilllabs.com

For HF27
leasing options or to
locate the Hill Dealer
nearest you,
please contact
Hill Laboratories and one
of our representatives will
personally assist you.