

# The Hill Hi-Lo



is designed to serve the doctor who uses the  
ACTIVATOR® instrument or applies the  
"ACTIVATOR METHODS®" technique.\* Featuring a smooth  
transition between soft start-up and landing,  
the Hill Hi-Lo will not disturb the patient with any sudden  
stopping or bumping. The standing foot board automatically  
folds down as the table top is lowered. Patient size and  
comfort are taken into account with a wedge-shaped foot  
section that conforms to the angle of  
the leg (also available with flat cushion).  
The Hill Hi-Lo is available in 20 decorator colors;  
the base is available in beige or gray.

## Specifications

Length: 6' • Width: 27"  
Optional Width: 30" (add \$50)  
Height: 34", optional 27"

**\$2995**



  
**Hill Laboratories**  
COMPANY

3 Bacton Hill Road, Box 2028, Frazer PA  
19355 610-644-2867 • fax 610-647-6297

[www.hilllabs.com](http://www.hilllabs.com)

\* (ACTIVATOR and ACTIVATOR METHODS are trademarks of Activator Methods, Inc. and are used for identification purposes only.  
Equipment made or distributed by Hill Laboratories is not made or approved by Activator Methods, Inc.)