

# PT Tilt Table

by Hill Laboratories  
2014

**The PT Tilt Table** gently moves your patient from a standing to a horizontal position. Featuring electric foot pedal-controlled movement and smooth transition between soft start-up and landing, the PT Tilt Table will not disturb the patient with sudden stopping or bumping. The PT Tilt Table is also equipped with three restraining straps that can be adjusted along the length of the table to fit different size patients. The PT Tilt is available in 20 decorator colors with base in beige, gray or black.

**\$3450**



## Standard Features

- Electric patient lift from standing to horizontal
  - Hospital Casters
  - 3 Restraining Straps
- Lifting Capacity - 400 lbs.
- Length: 6', Width: 27", Height: 34"
  - 2 1/2" medium-density foam
  - Choice of 20 colors for the top
- Beige or gray acrylic thermo-plastic base
- One-year warranty on all parts and labor

## Optional Features

- Width: 30" - \$50



The Hill PT Tilt Table moves your patient from vertical to horizontal position and back again smoothly and quickly.



**Hill Laboratories**  
COMPANY

3 Bacton Hill Road, P.O. Box 2028, Frazer, PA 19355  
Toll-Free 1-877-445-5020 • 610-644-2867 • fax 610-647-6297 • [www.HillLabs.com](http://www.HillLabs.com)